

# PRNC091A

Rhos

PRNC - 91-A

Revised June 18, 1970

fuente

PUERTO RICO NUCLEAR CENTER

MEDICAL SCIENCES AND RADIOBIOLOGY DIVISION

COURSE IN

TISSUE CULTURE AND RADIOISOTOPE TECHNIQUES

AT CELLULAR AND SUBCELLULAR LEVEL

PROGRAM AND LABORATORY EXERCISE MANUAL

June 18, 1970

1970

OPERATED BY UNIVERSITY OF PUERTO RICO UNDER CONTRACT

No. AT (W8-11-IEBD FOR U. & ATOMIC ENERGY COMMISSION

---Page Break---

ious

1

~PRNC - 91-0

COURSE IN

TISSUE CULTURE AND RADIOISOTOPE TECHNIQUES

AT CELLULAR AND SUBCELLULAR LEVEL

PROSRAM AND LAK

Bus

O\_o aR ese SS ee lo

---Page Break---

-\_? ee OS Oe oho ee oe oe oe: oe: ee ee ee ee ee oe ee oe i

---Page Break---

-ae

fur CURE eceRER,

Introdvetvien Lo Heme Cette Techniques

(Guaje 15 to 29, 1970)

Introductory Hem

re ~ Bistory of

je culture

oratory Bren

Peeparstion of Enterta),

Slerildzation and sleriity

Nests, Staining of Celis,

a

mp 16

\ 10:00 bats.

Feepuretien a

Constituents

weteeadey 17

itouredey 28

Heereine fs

Frahepors

of Cell Cultures

\_ Drs Je Chirivoga, Director

Medical Sciences and

Radiobiology Division

Dr. Re Hortfns

Drs Re Martfer-silva

Dr. Re Martfnensitve

Dre Re Hirtfaectivn

De. Ry Kertines

Dr. Re ar tinen-Stva,

De. Re itevtfnes- have

Dr. R, Mertinex-Siva

De. Re Martfaes Aiden

---Page Break---

a ae IS IS al ee. eee ae ae ee ee ae eee Se

Tuesday 23

Wednesday

Thursday 25

Friday 26

1s373 cuysuRs cou

une 15 to July 3, 1970

Vso of Teotoper at Celtuter and x

euler Te

Jocture labeling of Polysmes

Ye. Taboratory "OH

Tecture Muclate Aelia ant Protein

synthesis 1

|. Isborstory Inbo}ling of ?Polyscos (cont)

Tocture Metabolite yatlaays in cette



Use of Teotogee

10:00 git.

to 1:09 yn, Taberetry Metabolic

inaye An cos

Use of Teotepee

9:00 to scture Nucleic Acid end Protein sym

10:09 aan. thesis TZ

10:00 es, Inbora- Metabolic pathweya An cette

to 1:09 fa. tory, ?of Esotopes

9:00 to ecture ? Anatyice1 Wetheds 4a Cel

10:00 ya. end Holecular Biolay .

10:00 to Isboratory ??\_Centinuation

1:09 pa.

Dr.

Dr.

Dre

Dr.

Dre

Dre

Dr,

Dr.

Dr.

bie.

ana

2 (une 2 to 25, 1970)

Raymond Ae Bre:

Rayond Ay Brow

Jorge Chiriboga

Raynont A,

Tesvon-Aae

Roge

Roger Romo#-ALin

Foxge Chirivoe

Roper Raomh Lice

Raymond As Bs

Taymod A. Broun

Roger Hanoë-A4

---Page Break---

SS a eae eee rr eS Oe iim

lee eee ee ee ee

3

?Tissue CULTURE Cour

(une 29 to July 2, 1970)

Voulay 29 9:00 to Lgeture

39:00 am, he Belence of Virology

46

10:00 am

to a

1:00 pm, Measuring Virus Infection

by Cytopathic Effect

Laboratory

9:00 10 Lecture

10:00 am. Methods for Detecting Virus.

Multiplication in Time

culture

10:00 a.m. Laboratory exercise £9

Fleqie wettod for

Tayors of Celie

Lecture

Rethods for Detecting

Multiplication inv

culture

10:00 a.m. Inboratory Bueroine 3.19

to Woutralizing 4

2:00 pam, An Cold Cultures

?maveasy

to Lecture

10:00 a.m. ?Application of tesse Coitnre

in Virus Teolatie avd Voecine

Prouietion

Laboratory:

vation aud interpretstion

of Previous Fxporironts

lypLication of Sieeue Culture Techniques to Virology.

Dr, Julio I.

be. Julio Ts

be, dviso I.

Dr, dulso Te

De. Julio Te

be. Julio 1

Dr. Julio T.

Dr. Julio I.

coin

coro

conéa

coro

corn

corén

corn

coréo



---Page Break---

## LABORATORY PRACTISE #1

PREPARATION OF MATERIAL. STERILIZATION AND STERILITY

?TESTS. STAINING OF CELLS

aspirate

?ALL components in a cell or orga culture system must

be free of contaminant microorganisms, Routine testing

Procedures should be carried out in order to rule out the

Presence of bacteria, fungi and mycoplasma.

tests:

?two tubes of 1 culture

so these of test cells

Some tests of the 2 tests

Two tubes of tryptic soy broth

?Six tubes of Sabouraud agar .

five tubes of PLO bot

Solve tests of PO agar

Gram stain

Procedure:

Observe under the microscope and describe the different

cell types.

2 Tests for fertility

2:1 Inoculate 0.2 ml of each cell culture tube into 10 ml of

glycolate medium. Incubate at 37°C. Read and record

Results every 24 hours. If negative, discard after 5 days.

2:2 Streak @ loopfull of each cell culture on Sabourand Agar.

Incubate at room temperature, Read and record results

during 1 week, at 24 hour intervals,

---Page Break---

ach

Yahoratory Benetse #1

With @ pipette deposit some drops of each coll culture  
?eauta on Pio broth and agar (Bife0). streak the ager  
with @ bacteriological loop, invert the plate and  
Ancubate et 37° for seven days. Ineubste the broth  
?at the same texperature for four days and after this  
Period, place cose drops on a PPLO agar plate, spreading  
with a bacteriological oop.

Without removing the cover of the plates, inoculated  
?the 1st and lth days, Look for colonies of PPLO under  
?the microscope stage focusing through the agar. Use 4  
10K objective and 10, 12.5 or 15 x ocular.

Most FFLO colonies ce round, with  
?8 dense center end a less dense periphery, giving the  
appearance of a fried egg. YFLO colonies have been  
Aolated from tiseve cultures, however, they do not

conform strictly to this appearance on primary solution,  
?They may appear to lack a distinct periphery and appear  
to be totally embedded in the agar. ?these colonies are

usually very small and look "granular" or "fuzzy".

YELO colonies vary from 10 to 500 microns in diameter  
(0.02 to 0.5 mm) and characteristically the center only,  
Of all of the colonies, 4% embedded in the agar, Individual  
organisms cannot be resolved since they are the size of an  
average virus particle, Occasionally, at the periphery  
Of PELO colonies, "large bodies" characteristic of this

group of organisms are found, After locating the colonies,

---Page Break---

mBREeEe ee eS eS es as a2 eeUlcer EY OE eeOUDT,IIIOIU

Laboratory Exercise #1

they are usually marked out on the petri dish with a glass

?parking penest.

Confirmation of PELO colonies depends, in addition to

morphological characters, on:

1, Inability to remove the embedded portion of the colony from the agar surface by stroking the colony with a loop. This demonstrates the fact that the colony is embedded. Inert colonies will rub off.

2. The non-reversion to bacteria which subsequent passages of the colonies will reveal. Reversion to bacterial form would be typical of L forms.

3. A requirement for native protein.

4, Reaction with the Denes stain.

The Denes stain is prepared by dissolving 2.5 g

methylene blue, 1.25 g. azure II, 10.0 g maltose, and

0.25 g sodium carbonate in 100 ml of distilled water. With

cotton swab moistened in the stain, stroke the area of an

agar plate just adjacent to the suspected colony. The stain

will diffuse to the colony which is then examined under the

microscope as described above. The FLO colonies stand out

distinctly with densely blue staining centers and light blue

peripheries. Bacterial colonies are also stained but these  
are decolorized in about 30 minutes. The FLO colonies  
never decolorize the stain.

Staining of colonies by Gram's method.

Remove the growth medium of the tubes provided for sterility testing.

---Page Break---

Laboratory Exercise #1

2. Wash twice with Hank's solution.

3. Add methyl alcohol covering the surface of the slide. After 5

minutes the cellular sheet will be fixed.

Remove the alcohol and add 1 ml of Gram's stain (make a fresh  
stain by diluting 1 drop into 1 ml of distilled water).

5. After staining during 30 minutes, wash with tap water.

6. Remove the slide; let dry; mount following the instructions.

---Page Break---

8.

## LABORATORY EXERCISE #2

### PREPARATION OF TISSUE CULTURE CONSENTURNTS

Objective:

?To obtain the components which will provide the cell "in vitro"

with the organic and inorganic substances necessary for its nutrition,

and will protect it against changes in physical state, temperature,

pH, etc.

Materials: 7

4 chick embryo

Scissors

Syringe 50 ml.

Centrifuge tubes

Graduated cylinder

Antibiotics (Penicillin and streptomycin)

11, HANKS BALANCED SALT SOLUTION (R88)

A. 10 X Bolution

wait A Nalco; 3.5 am.

Discoive in 250 mi, diatilted water. Dispense in a convenient

Lottie (50 ml. screw-cep prescription bottle) and autoclave at

320°C for 15 minutes.

unit #2

wc 60.0 gn.

a 1.0 wn.

e904, 759 2.0 en.

a04.28,0 0.6 ox:

atucove 30.0 ex.

10370, 0.6 em

Pistolve in 800 ml, distitied vater.

---Page Break---



Laboratory Broth #2

unit #3

0.01 g.

Dissolve in 100 ml, distilled water.

unit #

Phenol Red O.M. am

Mix Phenol Red in a small amount of water until @ paste, dilute

to 150 ml. with distilled water, titrate to pH 7 with 1/20 acti.

Make up to final volume of 200 ml. Preserve with 1-2 ml, Chloroform.

MA 100 ml. of unit #2 and then add unit #3 to make

1,000 ml. Pour solution into glass stoppered bottle and add 3-4 ml.

chloroform as @ preservative. This solution may be kept at room

temperature for 6 months-1 year.

NOT: Minimize transfer of chloroform. In preparation of the working

solution. Be certain that bottle caps are loosened during

autoclaving to insure that all chloroform is driven off.

Working solution

The working BSS is prepared by diluting 10K Stock 1:10 with distilled water. Dispense in convenient size screw cap bottles and autoclave at 120°C for 15 minutes. Aseptically add 2.5 ml. of sterile sodium bicarbonate solution (Unit #1) to each 100 ml. of BSS. The pH may be adjusted with CO<sub>2</sub>. The balanced salt solution is now ready for use, Do not tighten cap until pH of BSS is 7.4.

## 2. NUTRIENT MEDIA

Eagle's minimum essential medium contains higher concentrations of amino acids than the basal medium first described by Eagle, which

---Page Break---

The medium may be prepared with Menck's PSS recipe.

Te mediua 18 prepared,

At the time of use,

red st 20°C) and NaHCO<sub>3</sub>; are added to the  
1X solution,

Solution Per Liter 10x negiun

L-Arginine. Hen

1.05 gn.

A-tystiaine. vo 0.33 em,

Atysine, Her 0.58 gm,

tryptophane 0.20 gm,

. 1-PhenyLatanine 0.32 an.

Lsthreonine © O48 gm,

Asteueine 0.52 gm,

I-Vatine 0.06 ee

I-Tsoleveine 0.52 em,

I-Methionine 0.15 om.

Solution B:

Antyrosine 0.32 gn.

aecyetine 0.24 gn,

?Tere anino acids are Aiseclved in 200 ml, of 0.075 Hea with  
gentle heating (60%),

Solution 0:

Meotinanide 200 me.

Pyridoxal, 200

?Thienine 200 ag

Pantothenic Acta 200 ng

Choline 200 ng,

SvInositor 400

Riboflavin 20 me.

---Page Break---

=

Laboratory Exercise £2

Components are aissolved in approxinetely:175 ml. of double  
aistinted water ana then brought to e final volume of 200 ml. with  
double distilled water. the solution 2 dispensed n 10 xi. ancunts  
and stored at -20°C. 10 ml, of Solution ø are added to each iter  
of 10x meatus.

solution Br

200 mi. of Biotin ere Gissolved in 150 ml. of double atstiiied

water. To increase stability during storage, 1.0 ml. of 1.0 HCL  
4s added. the total volume 4s brought to 200 ml. with double  
distilled water and the solution dispensed in 10 ml. aliquots and  
stored at -20°C. 10 ml. of Solution D are added to each Liter of  
20K medi.

Solution

200 mg. folie acta (crystalline) are dissolved in 200 ml. 1X Hanks!  
BES pH 7.8. The solution is dispensed in 20 ml. aliquots and stored  
at -20°C. 10 ml. of Solution E are added to each Liter of 10X medium,  
Glutamine Solution 36 - (To be added at the time of use) 12 ml. of

Asparagine are dissolved in 400 ml. of double distilled water and

sterilized by filtration through Seitz-type pad. The solution is  
stored at -20°C and 1.0 ml. is added to each 100 ml. of 1X Bagley's  
medium.

Preparation of the final mixture

iY BSS

fa, The following are dissolved in solution Bt

Wac 80.0 gx.

Ke 4.0 gee

2.0 gn.

---Page Break---

a2

Laboratory Exercise #2

B. The following are dissolved in 50 ml. double distilled water and

added to the pool.

apHPOy 12,0 1.52 em.

WLP, 0.60 on.

20 grams of Glucose are dissolved in 50 ml. of double distilled water with 20 ml. of 1% Phenol Red solution and added to the pool.

4. The volume of the pool is brought to 600 ml. with double distilled water and the following solutions are added:

Per 1.0 Liter 10x medium

Solution C 10.0 mL.

Solution D 10.0 ml.

Solution F 1 20.0 ml.

5. In 8 separate flasks containing 160 ml. double distilled water 2.0 grams of anhydrous  $\text{CaCl}_2$  are dissolved and added to the pool slowly with vigorous shaking.

6. The amino acids of Solution A are added to the pool and the volume is brought to approximately 950 ml. with double distilled water.

7. A solution containing 20,000 units of Penicillin in 20,000 micrograms of Streptomycin per ml. is added in a volume of 5.0 ml.

Per liter and the mixture is held in the refrigerator overnight.

The total volume is brought to exactly 2,000 ml. with double

distilled water and the solution is sterilized through a

type pad.

A. For use,

ø solution is diluted to 1X with sterile double distilled water and 24 of the 3f Glutamine Solution and 1.25 to 2.59, of @ 2.04 water, are added,

---Page Break---

## Laboratory Exercise #2

D. The following are dissolved in 50 ml. double distilled water and added to the pool.

NepHPO<sub>4</sub> 1290 1.52 on.

101%), 0.60 ex.

©. 20 grams of Glucose are dissolved in 50 ml. of double distilled water with 20 ml. of 1% Menol Red solution and added to the pool.

4. The volume of the pool is brought to 600 ml, with double distilled water and the following solutions are added:

Per 1,0 liter 10x medium

Solution ø 10.0 ml.

Solution D 10.0 ml.

Solution E 1 10.0 ml.



6. In a separate flask containing 160 ml. double distilled water 2.0 gms. anhydrous  $\text{CaCl}_2$  are dissolved and added to the pool slowly with vigorous shaking.

7. The amino acids of Solution A are added to the pool and the volume is brought to approximately 950 ml. with double distilled water.

8. A solution containing 20,000 units of Penicillin and 20,000 micrograms of Streptomycin per ml. is added in a volume of 5.0 ml. per liter and the mixture is held in the refrigerator overnight.

9. The total volume is brought to exactly 2,000 ml. with double distilled water and the solution is sterilized through a Seitz type filter.

10. For use, the solution is diluted to 1x with sterile double distilled water and 1% of the 3% Glutamine Solution and 1.25 to 2.5% of a 2.84  $\text{Na}_2\text{CO}_3$  are added.

---Page Break---

. \* laboratory Exercise #2

For Laboratories occasionally using smaller amounts of Page's media, it is recommended that the 1X medium be prepared by adding 4 ml 10X Stock Solution of the amino acids and the 20K stock solution

of the Vitenins (stored at -20°C) appropriately An Farle or Hanks?

299 adding Glutaaine, antibiotics and Nel003 as indicated shove.

?Mais prevents deterioration of the vitanine during? long term storage

at ice box temperature.

### 3. CELL DISPERSTNG AGENTS

A, Srypain Solution 1.0%

2 ge. of powdered trypein is dissolved in 100 ml. of phosphate buffer saline and the solution is passed through ash-free filter paper (Schleicher and Schull #589). ?The solution is then sterilized by filtration through a deite-type pad and stored at -20°C.

1B, Versene Solution (Bthyelenediamine tetraacetic acid)

Mech 8.0 om.

= HPO, 0.2 en.

Ke. 0.2 em

TiaghPO, 0.35 en.

Versene 0.20 gm.

Dissolve in 1,000 ail. of distilled water, Dispense in convenient

amounts and sterilize by autoclaving at 120°C for 15 mimites.

#### 4, AWETBIOTICS OLUTION

Fenieii2in (20,000 units per mi.) and Streptoryein (20,000 microgamta per mi.)

1. Ald 10 wd. of Marks" eolution to 1,000,000 unite of Penfci1iin,

2, Add 210 wl. of Honks' solution to 1 viel with a gram of Streptonycin.

---Page Break---

laboratory Frereise #2

3, Mix the contents of both vials end ada up to 50 ml. of Hanks! solution.

1, Dispense jo viels and keep at -20°C.

5. CHICK manyo EXIRACT 50%

Bach student will be provided with 5 exbryos 9-10 days of age.

A. Harvest the enbryos end place then in @ sterile Fetri-dish where the eyes, beake, lege and wings are removed.

B. The renaining tissvee are washed in a beaker containing Hanks! BSS, ?thon minced with uterine selssors.

©. the minced tiesve te passed through a 50 ml. ejringe into a

graduated cylinder or centrifuge tube.

D, An equal volume of Hanks' BSS is added to the tissue culture and the mixture is stirred and allowed to stand for 30 minutes.

B, the suspension is centrifuged at 1,500 rpm for 20 minutes and the supernatant fluid (constituting the 50% extract) is removed and stored at -20°C.

F. After thawing for use, the extract is clarified by centrifugation

at 2,000 rpm for 10 minutes.

---Page Break---

a5.

LABORATORY EXERCISE #3

MONOLAYER CULTURES OF CHICK EMBRYO CELLS

coast

satpind yates atures tren sick any were aoe

Sm frat to te nied for vine peopegation, evelomant of nor

testes hae Hoorensd the ute of chee eahryo firobst fr

fru fctaten ad engen geodetic.

Material,

i 4 enick enbryos

Petes dimes

Beaker

S| sefssors

Forceps

syringe

Erleaseyer flask

magnetic stirrer

Centrifuge tibes

Procedure

1. Chick exbryos 9 days old are harvested and placed in e sterile

Petri dish vere eyes, beaks, legs and wings are removed end aiscarded.

2: the embryos are transferred to beaker containing Ranke? BSS and

washed in 3 changes of the solution.

3. the extryos are minced Snto pieces approximately 3 mi. sn ataneter,

with uterine seiscore, and the minced tissue ie washed with 3

changes of Hanks? BSS.

4, Me minced tiasuc 4s passed through a 50 ml. syringe (without needle

into a 500 ml. Erlenmeyer Mask, vhere it 1e vashea twice with

90-100 ml. of Hanks' BSS.

---Page Break---

### Laboratory Exercise #3

5. After the fluid from the second washing has been removed, an appropriate volume of .25% solution of trypsin in Hanks' BSS is added. The volume of trypsin solution employed is determined by the number of embryos being processed (5-20 embryos 200 ml., 32-20 embryos, 300 ml.)

6. A magnetic stirring bar is added to the flask which is placed on the magnetic stirrer and the suspension agitated for 1 hour at room temperature,

After one hour of trypsinization the flask is slanted to sediment large tissue particles and the supernatant fluid (cell suspension) is decanted through a stainless steel wire cloth (72 mesh wire screen .0037 inches) into centrifuge tubes,

8. the cell suspension is centrifuged horizontally at 600 rpm for 20 minutes and the supernatant fluid containing trypsin is aspirated immediately. The cells are washed once by resuspending in Hanks' BSS, centrifuging at 600 rpm for 10 minutes and aspirating the

supernatant fluid.

9. The cells are then resuspended in 15 to 30 ml. of Hanks' BSS,

transferred to 15 ml. or 30 ml. graduated, conical centrifuge

tubes and centrifuged horizontally at 600 rps for 10 minutes.

10. The volume of packed cells is noted and after removal of the

supernatant

fluid, the cells are diluted in the following medium:

Povine serum 2.0 ml.

5 Lecithin hydrolyzate

in physiological saline 5.0 ml,

---Page Break---

an.

ae

Laboratory Exercise #3

50% Chick Embryo Extract 4.0ml.

Honks' Bss. 05.5 mi.

2.84 WalGO, 2.5 me

Penicillin-streptosyein solution 1.0 ml.

For tube culture the cells ere diluted 1 to 200 and dispensed in ml, volumes, For plaguing in stoppered bottles, the cells ere also diluted 1 to 200 end alepensed in @ xi. volumes into three ounces prescription bottles, For plaguing in 60 mm, Petri dishes 5 ml, of 1:200 ailution of the cells ere added and incubation is conducted in 9§ atmosphere.

After Sncubation at -37°C for 1 to 2 days, complete monolayers of celle are formed and the cultures sre ready for inoculation with viruses or clinical material.

---Page Break---

TADORATORY RKCESE gh

?TROHITLQUES FOR GROWING CELL LINES IN TISSUE CULTURE

crgectt

Tio experinent Ss designed to maintain a coll Line (tteta, ne 2 or 1)



for the duration of the course. Each student will be provided with 1 bottle of cells.

## Materials

1 bottle with a confluent cell monolayer

Haemocytometer

Solution trypan blue

14 tubes

Procedure: .

1. Observe the cells under the microscope and describe them.

2. Remove medium with a pipette. Wash

2 times with Hank's B99. Add

10 ml of 0.25 per cent trypsin solution, allowing the trypsin to cover the cells for exactly one minute at room temperature.

Remove all of the trypsin and place the tubes in 37°C

incubator for 10 minutes. At the end of this time the cells

should be almost completely detached from the wall of the bottle.

Add 10 ml of fresh medium and suspend the cells homogeneously

with a pipette,

3+ Count the cells in the haemocytometer

3-1, With a Pasteur pipette carefully express a drop of the cell

suspension made up of 0.5 ml. cells plus 1.0 ml. of trypan

Blue under the haemocytometer coverglass, avoiding any

overflow into the moat,

---Page Break---

laboratory Exercise 9

3-2, Determine the average number of viable cells (dead cells

\* stain (trypan blue) in the 4 large corner squares used for

counting white blood cells, multiply by 10,000 the

number of viable cells to obtain the number of cells

per ml. Adjust to 50,000 cells/ml. using Eagle's medium,

Transfer 50,000 cells into each of 1 test tubes. Stopper

tubes with rubber stoppers. Incubate at 37°C in an

horizontal plane.

4, Keep record and observe cells every any.

---Page Break---

## LABORATORY EXERCISE 45

### CONSERVATION AND TRANSFER OF CELL CULTURES

Objectives:

To maintain in the laboratory with a minimum of handling viable cell lines not in continuous use.

Materials

cell culture

Haemocytometer

Trypan blue solution

Aspores

A, Storage

1, Obtain a bottle with a culture of 5 day old cells from which the

out-growth medium has been removed

id replaced with 10 ml. of

fresh medium, consisting of 10 per cent horse serum and 90 per cent

Jactaltuain hydrolysate yeast extract medium. After a 2 day

Incubation period, the medium is removed, the cells trypsinized and

counted. Each bottle should yield between  $10^7$  and  $20 \times 10^8$  cells,

otherwise the cells are not suitable for storage.

2. The pH of the medium is adjusted to 7.4 by addition of an 8.0% HClO<sub>4</sub>

solution. Add 21.0 ml. of sterile glycerol to the 100 ml. of medium

4 ml each culture.

3. Transfer the cell suspension to ampoules (which can be freeze-dried)

or to tubes, tightening the stoppers and sealing with adhesive tape,

4. Bring to MC during 1 hour and then place the tubes (or ampoules) in

?the Reveo (~70°C) where the temperature will drop 1°C per minute,

Once the temperature inside the ampoule has reached a -20°C, the cells

---Page Break---

### Laboratory Exercise #5

can be stored at -197°C in the Liquid nitrogen refrigerator.

Under these conditions cells can be stored for periods up to 5 years.

5. To revive the frozen cells, the ampoule is removed from the Liquid nitrogen refrigerator and thawed rapidly in a 37°C water bath. A volume of the cell suspension containing  $1.5$  to  $2.0 \times 10^6$  is added to 10 ml. of outgrowth medium and cultures initiated in @ 200 ml. rotator.

#### 1B, Transport

1. Suspend a bottle of HeLa cells and dilute in growth medium to obtain a suspension of no more than  $0.6 \times 10^6$  cells/ml.

2, Refrigerate at 4°C for 24 hours,

3. centrifuge at 200 rpm for 30 minutes and discard the supernatant.

4. Add medium to obtain a cell suspension of  $2.4 \times 10^6$  per ml.

5. In this state the cells can be shipped in an iced container and

upon receipt sedimented by centrifugation at 200 rpm for 30 minutes

and resuspended in fresh growth medium at a concentration of  $0.6 \times 10^6$  cells/ml. A satisfactory method to ship cell cultures is obtained by placing the vessel with nutrient medium, preventing the trauma to the cells by the movement of the medium. Upon receipt, the medium should be removed and the cells fed with a new one.

---Page Break---

22.

## LABORATORY EXERCISE #6

### TUABELING OP POLYSONES

Objective:

To study the incorporation of RNA precursors into cellular polysones.

Material:

1 Liter RSB butter

0.02 M Tris-HCl pH 7.

0.01 M KCl

0.0015 M MgCl<sub>2</sub>

100 ml 10% sucrose w/w in SB

1100 ml 10% sucrose w/w in RSR

\ 20 ml ReB containing  $5 \times 10^{-7}$  g/ml hydrocortisone

500 ml Liquid scintillator .

2 flask bottles with confluent

monolayer of L-cells

? Trypsin for removal of cells

Procedure

Monday:

1. Make up isotonic solution.

2. Add 30  $\mu$ Ci uridine to each flask bottle.

3: take 6 10-204 sucrose gradients.

4, set up miniprep and flow cell for monitoring gradient

5. Set up paper discs for fractionation.

Tuesday:

1, Add 10  $\mu$ g/ml of uridine to Make bottles and incubate 10 min.

2, trypsinize cells,

---Page Break---

-23

### Laboratory Exercise #6

3. Centrifuge cells in conical centrifuge bottle.

4, Take up in 1.0 ml RSB containing  $5 \times 10^6$  g/ml hydrocortisone.

5. Homogenize in glass homogenizer.

6. Add 0.33 ml 20% sucrose to homogenate and 8000 RPM for 5 min.

7. Add 0.4 ml of BOC to supernatant and place 0.4 ml on each of 3 sucrose gradients.

8. Centrifuge in SW 39 rotor at 35,000 RPM for 30 min.



9. Pump gradients through Mow cell and on to proper aises.

20, Wash @iees 10 minut

in 10f TCA and 5 min. in 95% ethanot.

11, Iny @isce and adé to viele.

22. Count 1 min/viad.

---Page Break---

he

1 OF FOLYSOMES

Tis is a laboratory exercise, whose purpose is to fantltarize the student with cone techniques widely used in cellular and molecular bicology.

In the cell, proteins are synthesized on the polysomes. These consist of individual ribosomes held together by the messenger RNA which specifies the amino acid composition of the protein. The size of the polysome is proportional to the size of messenger RNA, consequently, of the peptide to be synthesized,

It would appear that the ribosomes are stable and can be reused by the cell whereas most messenger RNA is unstable. Hence the ribosomes are preferentially labeled by a long exposure to an RNA precursor (in this case [<sup>14</sup>C] uridine). Whereas the messenger RNA can be labeled by a short exposure to precursor ( [<sup>3</sup>H] uridine),

After breakage of the cell and removal of nuclei, mitochondria, and cell wall by low speed centrifugation, the principal particulate matter left in the supernatant is the polysomes. These are best analyzed by sucrose density gradient centrifugation. In this technique the sample is layered on top of a linear sucrose gradient in a centrifuge tube. There is a linear increase of sucrose with depth of the tube in order to stabilize the process of centrifugation.

After centrifugation in a swinging bucket rotor, it is necessary to analyze the contents of the tube. This can be done in various ways. However, we shall place the contents of the tube by injecting a more

dense liquid into the bottom of the tube. Total nucleic acid can be measured

---Page Break---

25

by allowing the displaced fluid to flow through a flow cell which is monitored at 260 mμ in a spectrophotometer. Fractions are collected on paper discs for analysis of radioactivity.

It is convenient to collect samples on discs since in this form they can be easily processed to yield defined biochemical fractions with uniform physical properties for counting (no necessity for quenching correction).

Liquid scintillation counting currently provides the simplest method for detection of β-emitters and the most frequently used technique for separation of double label. By an appropriate setting of the discriminators in to separate channels, both <sup>13</sup>C and <sup>14</sup>C can be counted accurately and efficiently.

---Page Break---

26.

?LABORATORY EXERCISE #7

METABOLIC STUDIES USING ISOTOPES IN TISSUE CULTURE

cagects

Uso of Innetnfonine-cMi, in the biosynthesis of phosphatiai? choline.

tncorporatfen of Cle for tranrdethiations mos

Materiade:

Hunan Liver cel cultures (Chang's cell Line)

Tamethionine-cUn, (Hew Mngland Tuclear Corp., Mecs.), sterile

Gatution (0-1 m2 2.5 ye).

?Teypein solution 0.25% in pit 7.4 phosphate buffer.

n-butencl. .

Teotenio sodium chloride solution (0.94).

Developing solvent for chronatograrhy:

+ ?ehloroform, methanol, water (65:25:

Gass plates with 250 uw Loyer of silicecgel G.

Preparation: Silicergei 6-15 ge," 30 ml distilled water; emlettied

and layered in each plate.

Chromatography chesbers.

Todine vaylor chanbers.

15 ml centrifuge tubes.

Sterile pipettes (0.5, 5.0 and 10.0 ml)

Scintillation Liquid (FPO 0.4f, FOROP 0.00f in toluene).

Radioactivity counting vials.

Beckman Liquid scintillator.

stirring rods and Pasteur pipettes

Proceai

1. Add 0.1 ml (2.5 µg) of L-methionine-<sup>35</sup>S dissolved in 10 µl of culture medium to a Chang's cell culture. Incubate for 2 hours

at 37°C.

---Page Break---

Laboratory Exercise #7

2. Trypsinize the cell culture with 5.0 ml of trypsin solution.

Resuspend the cells in 5.0 ml of culture medium,

3. Centrifuge at 1,000 RPM for 10 minutes.

4. Resuspend the cell pellet in 2 ml saline with gentle agitation.

Centrifuge at 1,000 RPM for 10 minutes,

Repeat the above procedure twice.

6. Resuspend the cells in 1.0 ml saline and transfer to a Total

Counting Vial with a Pasteur pipette, Rinse the  
counting tube with one addition of 0.5 ml saline portion,

7. Sonicate 3-5 minutes,

8. Add 1.0 ml of butanol. Shake for 2 hours,

9. Centrifuge at 3,000 RPM for 25 minutes.

20. Transfer 20% of the butanol phase to a counting Vial adding  
5.0 ml of the scintillating fluid.

UL. Prepare @ st24ca-get plate.

22. Dry in @ 200°C oven for 2 hours.

35. Allow the plate to cool and place an aliquot of the standard  
Phase 2 em from the Border, the volume depending on the  
Pedloactivity count, Next to this sample spot a standard of  
208 of Ine -tecithine,

2. Place 50 µl of the developing solvent in the chamber and then the  
stained plate,

15. Allow to develop until the solvent level has reached the top  
border.

---Page Break---

Laboratory Exercise #7

1G. Now to dry at room temperature, Develop the sample components  
in iodine vapor chamber,

1. Once the spots are identified the gel is collected from each one  
and placed in individual vials, Add 5 µl of scintillating liquid

found count in the Beckman Liquid scintillation counter using the carbon 14 tracer.

28. Discuss

of the results obtained.

---Page Break---

retrovirus a) Use of Stable isotopes in the study of cellular functions is based on the property of these elements to serve as molecular markers, whose fate is followed through the metabolic pathways. Stable isotopes are useful, the property of differentiation & state or radioactive source from the ordinary element, due to their presence in vast or by part of radiation, makes it

possible to determine the amount of the element in a sample.

Stable isotopes have the property of being incorporated into a marked compound within the cell or other level of higher organization, and their presence can be detected by the use of mass spectrometry.



guid scintillation system has a great advantage over others in the use of

radioactive emissions (# radiations of 13, etc),

D) Transmethylation mechanism to form phosphatidylcholine - this trans-

methylation process has been studied at different levels of organization.

of low energy

In intact animals, in homogenates, in cell fractions, in the presence of isolated enzyme systems, etc.

From intact animals (Du Vigneoul et al.), the first studies made on rats fed on a diet deficient in methionine and cysteine but in the presence of choline and homocysteine, led to a demonstration of the

existence of a metabolic process that forms choline from methionine.

---Page Break---

An evident proof of this operational mechanism was gotten "in vivo" using radioactive and stable isotopes. The use of  $^{14}\text{C}$ -methionine,

that, methionine and CDS-choline showed Ubet Were yas ay "ia toto!

Incorporating mechanism of the methyl groups transfer the methionine into the structure choline and viceversa and that the radioactivity was found in the phospholipide of the liver.

Other studies in liver homogenates and  $\text{S}_{100}$  fractions led to the identification of the same type of molecules involved in this process (Armer, J. and Greenberg, D.). In this way the radioactivity from  $^3\text{H}$ -methionine was incorporated in the phospholipide of these fractions and

$^3\text{H}$ -methionine in the presence of phosphoethanolamine, ATP, and  $\text{Mg}^{2+}$

cofactors formed phosphatidylcholine, phosphatidylethanolamine.

Previously, work dealing with  $\text{S}$ -Adenosylmethionine (Cantoni et al.) showed

methionine, ATP, and  $\text{Mg}^{2+}$  to be replaced by said compound in cell fractions

These works established the participation of at least two enzymes

that catalyze the formation of phosphatidylcholine by means of transmethylation

using methionine and a phosphatidic derivative. One that catalyzes the

a Me PF

(S-adenosyl methionine synthetase) and another that catalyzes the excessive

formation of S-adenosyl methionine using methionine, ATP,

methylation steps (S-Adenosyl methionine: phosphoethanolamine-ethyl transferase).

There is no information available on the study of this metabolic process in animal or human cell cultures, but the studies that have been made and

briefly published are valid for this level of organization, the cell level.

---Page Break---

CH, -0-0c-k'

© bH-e-0e-a"

e ø

BN = Cleo CI 0-P 0-8 Hy

ey ou Phrafatects whsvecle

aus Crene

whe CH ~CH~9 ~ pH ?O- CH,

Cle on

pie

BOW-ON CN. OP ~0-0 Hy

cu, by

Gus

i rt 1 ~O-

4 HCW, - Cl oP O-eu,

ey ow

---Page Break---

ae

## TARORATORY EXERCISE 48

Memions OF UETECTTHG OR MEASURING VIRUS MULTIPLICATION IN ISSU CULTURE

CHTOPATHIC EFFECT OF VIRUS MILITPLICATION

Qojectiv

4p determine the virus @i1ution that gives rise to cytopcthic

changes in 50 per cent of the inoculated eel cultures.

## Material

Poliovirus suspension

mbes

Pipettes

## Proceaur

<1, Observe the Heta cell cultures prepared during the third day of antoratory work.

2. Pipette off the grovth media aid replace St with .9 ml, each of Ragie's nedtun containing 2% serun.

3. Prepare tenfold attutions of polfovirus type 1 as follows:

a. Set up @ row of & Wassernon'e tubes munber

2 through ana

Aiepenee 1.8 ml, of media into each of then.

?b. Take 0.2 ml. of the poliovirus suspension end edd to the first,

?tube in the row. Mix thoroughly with a sterile pipette,  
c, Take 0.2 ml. from the dilution in the tube #1 and pass At to tube  
#2 mixing the contents.

@. Repeat the operation with the resainder tubes,  
4, mocutete .1 of esch dilution of the virus and deliver to each of  
4 tubes of lieta cells. A seperate pipette should be used for each  
Aiutions hovever in"the interest of laboratory glarsuare economy

use one 0.2 mi. pipette for adding the virus dilution to the Hela

---Page Break---

#### Laboratory Exercise #8

cette beginning with the highest Afiution and working back to the  
dover.

5. Set up to 2 Fela tubes for control without inoculation.
6. Bring to the incubator at 37°C and read record the results every day.
- T. Calculate the TEIgy by the Reed-ttiench method.

In the following table an exanple ir given of data derived from an

ideal experizent for illustrating the procedure of secumlation:

we

7 ?ACCUMULATED VALUES

a a cs

dilution ' Ratio CPE ?Wo. CFE \_' CPE "No. CPs" Ratio "Per Cent

Wife tw "0 'a2 ' 6 ' 32/2 + 100

wes Wr w+ 6 § 8 + ot BB + 20

ws a gt a kh at st oo

ao ahh. \* 3B fat kt apt 20

wF + of + o \* rr ee °

Accumated velues for the total muaber of tybes that shoved a CPB  
or were intact are obtained Ly adding in the direction indiested by the  
?srrows, the accumuliated CTE ratio represeste the acoumiated maber of  
tubes with eytopathhte changes over the secumated total mber inoculated.  
Jn this example the cytopathic change in the 10°3 asiution, 18  
higher than 50%5 4m tho next lower @itution, 1074, it is considerably

Jower. the necessary proportionate distance of the 50 per cent CPE end

---Page Break---

-33+

Taboratory Brereise £8

8 obtained as follows:

oint Lies between these tvo dilutions end

Per cent CPE at dilution next above 50f) ~ Sof = Proportionate

Per cont CPE at dilution next above 50%) - (per cent CPR \_adatanca

?at ablution

next below)

or 80 - 50: 30 x 0,5

B20 bo

Since logarithmically the distance between two ailutions is @

funtion of the incremental steps used in preparing the series, st is

necetsary to correct the proportionate aistance by the dilution factor.

In the case of seriat ten-fold dilution the factor 4s 1 (10g 10 © 2) and

Yuin is disregarded. In our exemple we have:

Negative logarithm of CTD, and point titer = negative logarithm of the

D,



30 dilution above the 50 per cent CPE

\ plus the proportionate distance

\ F factor

that is:

Negative logarithm of the dilution above 50% -3

Proportionate distance (0.5 X dilution factor (log 10): 0.5

TCTs 35

= 10<sup>-35</sup>

Log TCM = 10

---Page Break---

he

LABORATORY EXERCISE 49

PLAQUE METHOD FOR CONFLUENT LAYERS OF CELLS

Objective:

To produce circumscribed infected areas by vaccinia virus in chick

embryo fibroblasts which do not take the neutral red vital stain

some appear as clear unstained areas against a background of viable stained cells,

waterbath

Calf embryo fibroblasts

Vaccinia virus suspension

Water bath

Neutral reaction

Petri dishes :

Procedure

1. Set up three Wasserman's tubes numbered 1 through 3 and add 1.8

ml. of Eagle's media to all tubes. From the pool of vaccinia virus

supplied add .2 ml. to tube #1

Mix thoroughly with sterile

Pipette, Withdraw 0.2 ml. and add to tube #2. Repeat the previous

step with new pipettes so that you will have virus titrations 10<sup>7.2</sup>,

10<sup>7.3</sup>, and 10<sup>7.4</sup>,

3. Remove the out-growth medium from the Petri dishes previously Prepared with chick exbryo fibroblasts and vush once with sxline. Inoculate a plate with 0.5 ml. of each virue dilution.

3. Incubate at 37°C for 2 houre for virve adeorption with occastonal

rocking to distribute the virus particle:

4. After the adsorption period, remove the fluid ?and overlay the fibroblast sheets with 5 mL. of the following agar medium:

Agar 1.0 gm.

---Page Break---

1

Yeast Bxtract (Difco) oa em

actalbunin hydrotysate oh am.

Horse seri 34.0 mi.

Take? Bee 85.0

Fenici22in 50 u/s. end 50 Streptomycin to complete medium.

5. Allow the agar to solidify and tum the plate upside down and

Ancubate at 37°C for 3h days with the cell monolayer dow.

6. Add 3 ml. of a 1/1000 1+ 7 neutral-red solution; incubate at room

?temperature for 2-h heure and overnight at °C,

7. The dilutions used should produce distinct end separated plaques.

Observe against a vbite background. Hy counting the munber of

plaques at the Ailution sere they appear distinct, and by

multiplying by the cocr-cpondent ?idution fector, the number of

plaque forming units (PF) por ml. of the virus suspension can

be ealeuisted.

---Page Break---

TABORATORY EXERCISE 410

UBUYRALT2ING ANDTBODIRS ASSAYS IN CELL CULTURES

Sajectiver

to eteviate the cepeeity of a serum to neutralize the eytopathse

effect of poliovirus in a HeLa cells systems.

Materials:

Serum

Poliovirus suspension

Subes with HeLa ces

Water bath at 56°C.

Procedure:

1. the serum specimen is inactivated at 56°C for 30 minutes to

destroy heat non-specific virus inhibitory substances.

2. Set up a row of tubes to make serum dilutions of 1:4, 1:16,

1:64, 1:256,

1:1024 prepared in either balanced salt solutions

or the maintenance medium to be used in the cell cultures,

3. Poliovirus diluted to contain 100 TCID<sub>50</sub> & volume of

0.1 ml, (as determined by a previous titration of the virus).

The viral dilutions are made in the same medium employed for

the preparation of the serum dilutions.

4, Equal volumes of the serum dilutions (0.5 ml.) and of the diluted

test virus (0.5 ml.) are mixed, the volume of serum virus mixture

prepared is dependent upon the number of cell cultures to be

inoculated with the mixture. For virus control the test virus

dilution is mixed with an equal volume of diluent (or known

normal serum) and incubated under the same conditions as the serum

---Page Break---

For laboratory reference #20

virus mixture. For serum control, the 1/4 dilution is mixed

with antigen, it is necessary to perform a concurrent titration

of the virus to establish that a test dose actually contains approximately 100 TTB.

The conditions recommended for incubation of serum virus mixtures

vary widely for certain agents. It has been demonstrated that some preliminary incubation does increase the neutralizing capacity of serum. The important consideration is to avoid incubation conditions under which the virus may be labelled for a sample for long periods at 37°C. For most neutralization tests; incubation of the serum virus mixtures is conducted for 30 minutes to one hour at room temperature or at 4°C.

[After incubation period, serum virus mixtures, virus controls and serum controls, are inoculated in aliquots of .2 ml. into monolayer tube cultures: At least two cultures are employed for each serum mixture «

the inoculated cultures are incubated at 37°C and examined microscopically for ability of the serum to inhibit CPE of the virus.

The cytopathic effect of the inoculated tubes is recorded

find the neutralizing end point ie expressed av that elution

of serum which protects 50% of the tubes against the test doce

---Page Break---

ER EEE EE SE'SSE\_? eo eel! IL

38

laboratory Exeroise #10

of the virus, as {lustedated in the following examie:

sem DIRION \* CRE? CRE "Ho cHPE + ope ino con | MORTALITY

ae ?sotto 1 1 MOS patio ter cont

uh Foe to 2 te 6 toe °



me ige fe la le la im los

a:6h tye tas a ta te tap B

wee tue fa fa ba fa tap lg

10h app eo tk te Am + 100

SSPE gee oy os

logarithm 50 per cent stralizing end point =  $-1.8 + (0.5 y (0.6)$

«28+ 60.)

= 2a

antilogaritim -2.1 = 120

---Page Break---